

## Zebrafish spinal canal cilia staining

Mia Konjikusic

- Euthanize embryos in high-dose Tricaine and fix in 4% sweet PFA (4% PFA, 4% sucrose in PBS) overnight at 4 °C.
- Wash embryos in PBSTr (1x PBS with 0.1% Triton X-100), 3x and distilled water 2x 5 min. each
- Incubate embryos in acetone for 7 min at -20 °C
- Wash again in distilled water 2x and in PBSTr 3x 5 min each
- Incubate in PBD blocking solution (1% BSA, 2% DMSO, 2% normal goat serum in PBST) for 1 hour
- Incubate embryos with a mouse GT335 glutamylated tubulin antibody primary antibody (1:500) diluted with PBD overnight at 4 °C
- Wash in PBD for 2x 10 minutes each then washed in PBSTr 3x 10 minutes each
- Add secondary antibodies (1:1000, diluted with PBD) overnight at 4 °C.
- Wash in PBD for 2x 10 minutes each then washed in PBSTr 5x 10 minutes each
- Mount embryos in low melt agar and image

*Never remove all the liquid during washes, keep embryos fully submerged.*

*Washes and antibody staining can go longer but extended times can lead to deterioration of the tissues.*