

Colorimetric and Florescent In Situ Hybridization

Modified from Mariani lab

Zhaoyang Liu 2017.6

Day1

1. Deparaffin slides

Xylene	5min
Xylene	10min
100% EtOH	5min
95% EtOH	5min
70%EtOH	5min
1XPBST (0.1%tween-20)	5min

2. Permeabilization

5min incubation with 7.5ug/ml proteinase K (in PBST) (on slide)
5min wash in PBST in a Jar (1)
Postfix in 4% PFA with 0.2% Glutaraldehyde for 20min (on slide)
(add 8ul 25% glutaraldehyde in 1ml 4% PFA)
5min wash in PBST Jar (1)

3. Hybridization
 - a. Pre-hybridise slides in **preheated** hybridization buffer (HYB buffer) at 68C for (at least) 1h in a Jar (2).
 - b. Dilute the RNA probe to 10ng-20ng/ul with HYB buffer, heat on the heat board at 80C for 10min, and then put directly on ice. Apply 100ul probe per slide. Cover with a cover slide and then cover with a piece of para film.
 - c. Put slide in a humidified chamber with filter paper soaked with H2O or HYB buffer at the bottom. Wrap the chamber if necessary. Hybridise with probes overnight (12-16 hours).

HYB buffer (200ml):
100ml Formamide, 50ml 20XSSC (ph 7.0), 1g Torula Yeast (tRNA) (**R6625 SIGMA**), 10mg Heparin Sodium Salt (50ug/ml), 50ml H2O.

20X SSC (500ml):
87.65g Sodium Chloride, 44.1g Sodium Citrate, dissolved in H2O and adjust the pH to 7.0 with Citric acid and bring the solution to the final volume of 500ml.

Day 2 (Florescent)

4. Stringency washes
 - a. Wash slides in **preheated** wash buffer at 68 1X15min, 2X30min in a jar (2)
Do not remove the cover-slide from the slide. Place the slide with the cover-slide into the jar and let the cover0slide detach from the slide in the buffer.
Wash buffer (200ml)
100ml Formamide, 10ml 20XSSC, 200ul Tween-20, 89.8ml H2O.
 - b. Wash slide in TNT buffer 3X15min at room temperature in a Jar (3).
TNT (1 litter): 100ml 1M Tris pH7.5, 30ml 5M NaCl, 1ml 100% Tween-20, 869ml H2O.
5. Blocking endogenous Peroxidase
Incubate in 2% H2O2 diluted in TNT for 30-60min at RT in a Jar (3).
Add 2ml of 30%H2O2 to 30ml TNT.
Wash in TNT for 1X5min in a Jar (3).
6. Blocking in 2% Roche Blocking Reagent (11096176001 Roche, in 1XTNT) for at least 1h on slide at RT.
7. Antibody incubation.
Dilute antibody with 2% Blocking Reagent: (start with the strongest probe)
Anti-DIG-POD 1/500 (11207733910 ROCHE, Anti-Digoxigenin-POD, Fab fragments)
Anti-DNP-POD 1/200 (Perkin-Elmer Cat# NEL747A001KT)
Anti-Fluo-POD 1/200 (Invitrogen cat# A-21253)
Incubate at 4C overnight (on slide). Cover the section with a piece of para film.

Day 3 (Florescent)

8. Color reaction
 - a. Wash 3x5 min in TNT in a Jar (3)
 - b. Equilibrate sections for 5min with Perkin Elmer Amplification Diluent (supplied with the kit) (on slide) (TSA Cyanine 3 and Fluorescein system: Perkin Elmer Cat# NEL753001KT)
 - c. Prepare appropriate dilution of the fluorophore in amplification diluent. (1:100 in amplification diluent)
 - d. (the following steps should be prepared in dark to avoid bleaching)
Apply 100ul of antibody per slide; incubate for no more than 1h at RT. Cover the section with a piece of para film.
 - e. Wash in TNT 1x5min in a Jar (3)
 - f. Quench reaction by washing in 2% H2O2 in TNT for 30-60min. (If this step is incomplete, you would get cross-reaction with the following staining. Can skip if you only have one color)
 - g. Wash in TNT 1x5min in a Jar (3).
 - h. For developing the next color, repeat step 6-8 .
9. Nuclear counterstain: counterstain with DPAI (300nM) for 5min at RT, wash with 1X5min PBS in a Jar (4).

10. Mounting with mounting medium.

Day 2 (Colorimetric)

4. Stringency washes
 - a. Wash slides in **preheated** wash buffer at 68 1X15min, 2X30min in a jar (2)
Do not remove the cover-slide from the slide. Place the slide with the cover-slide into the jar and let the cover0slide detach from the slide in the buffer.
Wash buffer (200ml)
100ml Formamide, 10ml 20XSSC, 200ul Tween-20, 89.8ml H2O.
 - b. Wash slide 3X15min in MABT in a jar (3)
MABT (1litter): 200ml 5XMAB, 1ml 100% Tween-20, 799ml H2O
5X MAB: 500mM Maleic Acid, 150mM NaCl, adjust to pH7.5 with NaOH.
5. Blocking in 2% Roche Blocking Reagent (11096176001 Roche) for at least 1h on slide at RT.
6. Antibody incubation.
Dilute antibody with 2% Blocking Reagent
Anti-DIG-AP 1/2000 (11093274910 ROCHE, Anti-Digoxigenin-AP, Fab fragments)
Incubate at 4C overnight (on slide). Cover the section with a piece of para film.

Day 3 (Colorimetric)

1. Color reaction
 - a. Wash 3x5 min in MABT in a Jar (3)
 - b. Wash in staining buffer NTMT 2X5min in a Jar (4)
NTMT(100ml): 10ml 1M Tris pH9.5, 5ml 1M MgCl2, 2ml 5M NaCl, 100ul 100% Tween-20.
82.9ml H2O. The Color reaction is highly dependent on the pH. Make sure the pH of the buffer is 9.5.
 - c. Apply staining solution **BM purple** to each slide. Cover the section with a piece of para film.
Keep the slide in dark for 2h-overnight at RT. After one day, you can change the staining solution. Wash the slide in NTMT and they apply fresh staining solution.
 - d. Wash 3X5min in PBST.
2. Counterstaining
15 dips in H2O
1min incubation in nuclear fast red
15 dips in H2O
10min incubation in H2O in a Jar (4)
3. Dehydrate with 2X2min 95% EtOH
2X2min 100 EtOH
2X2min Xylene
Mounting with Cytoseal 60 and cover slide.